

```
//#include "FitApeak.hh"
#include "Fit2peaks.hh"
#include "TCanvas.h"
#include "TH1D.h"
#include "TFile.h"
#include "TTree.h"
#include "TString.h"
#include "TAxis.h"
#include "TStyle.h"
#include "TSystem.h"
#include <fstream>
#include <utility>

#include <iostream>
using std::cout;
using std::endl;

#include <vector>
using std::vector;

#include <string>
using std::string;

#include <sstream>
using std::stringstream;

int main(int argc, char* argv[])
{
    int runNumber;
    string inputCalName("");
    string outputCalName("");
    string reportName("");
    string resultsName("../output/mg25Results%i.txt");

    if(argc==1)
    {
        runNumber = 244;
        inputCalName = "../cal/calDefault.txt";
        outputCalName = "../cal/new/calDefault.txt";
        reportName = "../reports/test.pdf";
    }
    else if(argc==2)
    {
        stringstream ssRunNumber(argv[1]);
        ssRunNumber >> runNumber;
        inputCalName = Form("../cal/cal%i.txt", runNumber);
        outputCalName = Form("../cal/new/cal%i.txt", runNumber);
        reportName = Form("../reports/report%i.pdf", runNumber);
    }
    else if(argc==3)
    {
        stringstream ssRunNumber(argv[1]);
        ssRunNumber >> runNumber;
        inputCalName = Form("%s", argv[2]);
        outputCalName = Form("../cal/new/cal%i.txt", runNumber);
        reportName = Form("../reports/report%i.pdf", runNumber);
    }

    // remove annoying "pdf/png created statements"
    gROOT->ProcessLine("gErrorIgnoreLevel = 1001;");
    gStyle->SetOptStat("ni");
}
```

```

gStyle->SetOptFit(1111);

//run10427.root is background run at the end of the low energy portion of the
2017 170(a,ng) HAGRID runs
//
const int Ndetectors = 13;
const int bnum[Ndetectors] = {0,0,0,0,0,0,0, 1,1, 1,1,1,1 };//
const int cnum[Ndetectors] = {0,1,2,3,4,5,6, 0,1, 3,4,5,6 };//

TCanvas *c1 = new TCanvas("c1","c1",0,0,600,400);

//runfile name pattern for COMPASS
TString runfilename = Form("../rootfiles/run_%i/UNFILTERED/
compass_run_%i.root",runNumber,runNumber);

cout << "Trying to find UNFILTERED Data... " << runfilename << endl;
if(gSystem->AccessPathName(runfilename.Data() ) )
{
    cout << "\tUNFILTERED Data not found!... " << runfilename << endl;
    runfilename = Form("../rootfiles/run_%i/FILTERED/
compass_run_%i.root",runNumber,runNumber);
    cout << "Trying to find FILTERED Data... " << runfilename << endl;
    if(gSystem->AccessPathName(runfilename.Data() ) )
    {
        cout << "\tFILTERED Data not found!... " << runfilename <<
endl;
        return -1;
    }
}
cout << "\t Processing data file " << runfilename << endl;
TFile *runFile = new TFile(runfilename);
TTree *runTree = static_cast<TTree*>(runFile->Get("Data"));
TFile *bgfile = new TFile("../cal/calibratedBackgroundHists.root");
Fit2peaks *mg25analysis = new Fit2peaks(bgfile);
mg25analysis->SetEgamma(0,1778.97);
mg25analysis->SetEgamma(1,2838.3);
c1->Print(Form("%s(",reportName.c_str()),".pdf");
ifstream ecal(inputCalName.c_str());
ofstream ecalNew(outputCalName.c_str());
ofstream results[Ndetectors];

//mg25analysis->Rebin(16);
for(int i = 0; i<Ndetectors;i++)
{
    results[i].open(Form(resultsName.c_str(),i),std::ofstream::app);
    if(!results[i].is_open())
    {
        cout << "error opening results file: " <<
Form(resultsName.c_str(),i) << endl;
    }
    if(!ecal.is_open())
    {
        cout << "error opening calibration file: "
<< inputCalName << endl;
        break;
    }
    TString id,check;
    check = Form("b%i_c%i",bnum[i],cnum[i]);
    double offset,gain,quad,resolution;
    ecal >> id >> offset >> gain >> quad >> resolution;

```

```

        if(id!=check)
        {
            cout << "channel mismatch with caibration (" +id+"!="+check+")
"
            << inputCalName << endl;
            break;
        }
        //T0 D0: grab this time from run.info file
        double tfg = runTree->GetMaximum("Timestamp")-runTree-
>GetMinimum("Timestamp");
        tfg/=1e12; // Timestamps are in picoseconds.
        double tbg = 3987.0;// Background run live time in seconds from the
timestamps in run 244

        mg25analysis->LoadBackground(bnum[i],cnum[i]);
        mg25analysis-
>CreateSpectrum(runTree,bnum[i],cnum[i],Form("h%i_%i",bnum[i],cnum[i]));

        mg25analysis->SetParameters(1000,resolution,offset,gain,quad,
1,0.002,500,1000,resolution*1.5,0.002);
        mg25analysis->FixParameter(5,tfg/tbg);
        c1->cd();
        c1->SetGridy();
        c1->SetLogy(1);
        mg25analysis->Reject(true);
        mg25analysis->GetSpectrum()->Draw("E0");
        mg25analysis->GetSpectrum()->GetXaxis()->SetRangeUser(0,0);
        mg25analysis->Fit("mq");
        mg25analysis->Subtract();
        mg25analysis->Reject(false);
        mg25analysis->GetBgFunction()->DrawClone("SAME");
        mg25analysis->Reject(true);

        //mg25analysis->RejectMore(true);
        mg25analysis->GetSpectrum()->GetXaxis()->SetRangeUser(250.0/gain,
5000.0/gain);
        c1->Update();
        c1->Print(reportName.c_str(),".pdf");
        mg25analysis->GetSpectrum()->GetXaxis()->SetRangeUser(1400.0/gain,
2200.0/gain);
        c1->Update();
        c1->Print(reportName.c_str(),".pdf");
        mg25analysis->GetSpectrum()->GetXaxis()->SetRangeUser(2500.0/gain,
3500.0/gain);
        c1->Update();
        c1->Print(reportName.c_str(),".pdf");
        mg25analysis->GetCalibrated()->Draw("E0");
        mg25analysis->GetCalibrated()->GetXaxis()->SetRangeUser(400,3000);
        c1->SetLogy(0);
        c1->Update();
        c1->Print(reportName.c_str(),"pdf");

        mg25analysis->GetCalibrated()->GetXaxis()-
>SetRangeUser(1779-150,1779+150);
        c1->Update();
        c1->Print(reportName.c_str(),"pdf");
        gain = mg25analysis->GetParameter(3);
        offset = mg25analysis->GetParameter(2);
        resolution = mg25analysis->GetParameter(1);
        double yield = mg25analysis->GetPeakArea(0);
        double error = mg25analysis->GetPeakUnc(0);

```

```
        ecalNew << id << "\t" << offset << "\t" << gain << "\t" << quad << "\t"
<< resolution << endl;
        results[i] << runNumber << "\t" << yield << "\t" << error << "\t" <<
gain
        << "\t" << offset << "\t" << resolution << endl;

    }
    ecal.close();
    c1->Clear();
    c1->Print(Form("%s)", reportName.c_str()), ".pdf");

    delete mg25analysis;
    delete bgfile;
    delete runFile;
    return 0;
}
```